

The kit **Enzy-Well Treponema IgG (REF 91050)** is an immunoenzymatic method for the qualitative determination of IgG-class antibodies to Treponema pallidum in human serum.

Each kit is composed of:

<b>MT PLATE</b>	<i>MICROPLATE</i> <span style="float: right;"><i>12 x 8 wells coated with a pool of recombinant proteins of Treponema pallidum</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet.
<b>CONTROL +</b>	<i>POSITIVE CONTROL</i> <span style="float: right;"><i>1 x 1.6 mL</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet. <u>Contents:</u> diluted human serum, containing anti-Treponema pallidum IgG antibodies, in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%.
<b>CONTROL CUT OFF</b>	<i>CUT-OFF CONTROL</i> <span style="float: right;"><i>1 x 2 mL</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet. <u>Contents:</u> diluted human serum, containing anti-Treponema pallidum IgG antibodies, in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%.
<b>CONJ</b>	<i>CONJUGATE</i> <span style="float: right;"><i>1 x 16 mL</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet. <u>Contents:</u> anti-human IgG monoclonal antibodies labeled with peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%.
<b>CONTROL IgG -</b>	<i>IgG NEGATIVE CONTROL</i> <span style="float: right;"><i>1 x 1.6 mL</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet. <u>Contents:</u> human serum not diluted in phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%.
<b>WASH BUF 10x</b>	<i>WASHING BUFFER 10X</i> <span style="float: right;"><i>1 x 100 mL</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet. <u>Contents:</u> Phosphate buffered saline (PBS), concentrated 10 times; contains Brij 0.5%
<b>SAMP DIL 10</b>	<i>DILUENT 10</i> <span style="float: right;"><i>1 x 100 mL</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet. <u>Contents:</u> Phosphate buffered saline (PBS), containing proteins 10%, sodium azide 0.09% and methyl orange.
<b>SUBS TMB</b>	<i>SUBSTRATE</i> <span style="float: right;"><i>1 x 12 mL</i></span> <u>The substance is not classified as dangerous according to Regulation 1272/2008 and does not contain hazardous substances and/or with limit values for exposure</u> in concentrations above the established limits. For this mixture, according to EC Regulation 1907/2006, is not necessary to develop a material safety data sheet. <u>Contents:</u> Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8).
<b>H<sub>2</sub>SO<sub>4</sub> 0.3 M</b>	<i>STOP SOLUTION</i> <span style="float: right;"><i>1 x 16 mL</i></span> <u>The substance is classified as <b>dangerous</b> according to Regulation 1272/2008.</u> <u>The material safety data sheet, prepared in accordance with EC Regulation 1907/2006, is available below.</u> <u>Contents:</u> solution of sulphuric acid, 0.3 mol/L.



**MATERIAL SAFETY DATASHEET**  
**ENZY-WELL STOP SOLUTION**  
(According to Regulation (EC) No. 1907/2006)

Edition: 4  
Date: 18-Jun-2018  
Page: 2/7

**1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND THE COMPANY/UNDERTAKING**

**1.1 PRODUCT IDENTIFIER**

**Product name:** ENZY-WELL STOP SOLUTION  
**Product number:** PF 93602  
**Brand:** DIESSE

**1.2 RELEVANT IDENTIFIED USES OF THE SUBSTANCE/MIXTURE AND USES ADVISED AGAINST**

**Identified use:** Professional use as laboratory reagent.

**1.3 DETAILS OF THE SUPPLIER OF THE SAFETY DATA SHEET**

**Company:** DIESSE Diagnostica Senese SpA

Registered office:  
Via A. Solari 19  
20144 Milan, Italy  
Tel: +39 02 4859121  
Fax: +39 02 48008530

Production plant:  
Via delle Rose, 10 – Loc. Tognazza  
53035 Monteriggioni (SI), Italy  
Tel: +39 0577 587111  
Fax: +39 0577 318690

**1.4 EMERGENCY TELEPHONE NUMBER**

**Emergency number:** Centro Antiveleni, Ospedale Niguarda Ca' Granda - Milano  
Tel: +39 02 66101029

Centro Antiveleni, Azienda Ospedaliera "S.G.Battista" – Molinette di Torino - Torino  
Tel: +39 011 6637637

Centro Antiveleni – U.O. Tossicologia Medica, Azienda Ospedaliera Careggi – Firenze  
Tel: +39 055 4277238

Centro Antiveleni, Policlinico A. Gemelli – Università Cattolica del Sacro Cuore - Roma  
Tel: +39 06 3054343

Centro Antiveleni, Azienda Ospedaliera A. Cardarelli – Napoli  
Tel: +39 081 7472870

web: [http://www.salute.gov.it/servizio/documenti/centri\\_antiveleni.pdf](http://www.salute.gov.it/servizio/documenti/centri_antiveleni.pdf)

**2. HAZARDS IDENTIFICATION**

**2.1 CLASSIFICATION OF THE SUBSTANCE OR MIXTURE**

**Classification according to Regulation (EC) No. 1272/2008:**  
Skin Corrosion – Category 1A

**Dangerous components:**  
H314 - Causes severe skin burns and eye damage

**2.2 LABEL ELEMENTS****Pictograms:** GHS05**Signal word:** Danger**Hazard statement:**

H314 - Causes severe skin burns and eye damage

**Precautionary statement:***Prevention*

P260 - Do not breathe dust/fume/gas/mist/vapours/spray

P264 - Wash ... thoroughly after handling

P280 - Wear protective gloves/protective clothing/eye protection/face protection

*Reaction*

P301 + P330 + P331 - IF SWALLOWED: rinse mouth. Do NOT induce vomiting

P303 + P361 + P353 - IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower

P363 - Wash contaminated clothing before reuse

P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

P310 - Immediately call a POISON CENTER or doctor/physician

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

*Storage*

P405 - Store locked up.

**Contains:**

Sulphuric Acid

INDEX: 016-020-00-8

**2.3 OTHER HAZARDS**

None

**3. COMPOSITION/INFORMATION ON INGREDIENTS****3.1 SUBSTANCE** Non applicable**3.2 MIXTURE**

Substance name	Concentration	EC Regulation 1272/2008
Sulphuric Acid		
Cas No 7664-93-9	<3%	Skin Corr. 1A H314
EC No 231-639-5		
Index No 016-020-00-8		

The entire text of Hazard statements is reported ay section 16 of the sheet.

**4. FIRST AID MEASURES****4.1 DESCRIPTION OF FIRST AID MEASURES****General advice:**

Move the person from the exposure to open air. In case of severe condition consult a physician immediately and present this safety data sheet. Eyewash and shower for the treatment of emergency have to be present in the workplace.

**Skin contact:**

Take off the contaminated clothes and wash with copious amounts of water. If irritation persist consult a physician.

**Eye contact:**

Take off lenses (if present). Flush eyes with plenty of water for at least 15 minutes keeping eyelids opened. Consult immediately a physician.

**Inhalation:**

In case of severe condition consult a physician.

**Ingestion:**

Clean mouth with water and drink afterwards plenty of water. Do not induce vomiting. Do not administer any medicine if the person is unconscious. Consult immediately a physician.

**4.2 MOST IMPORTANT SYMPTOMS AND EFFECTS, BOTH ACUTE AND DELAYED**

No data available



**MATERIAL SAFETY DATASHEET**  
**ENZY-WELL STOP SOLUTION**  
(According to Regulation (EC) No. 1907/2006)

Edition: 4  
Date: 18-Jun-2018  
Page: 4/7

**4.3 INDICATION OF ANY IMMEDIATE MEDICAL ATTENTION AND SPECIAL TREATMENT NEEDED**

No data available

**5. FIREFIGHTING MEASURES**

**5.1 EXTINGUISHING MEDIA**

**Suitable extinguishing media:** Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

**Unsuitable extinguishing media:** None

**5.2 SPECIAL HAZARDS ARISING FROM THE SUBSTANCE OR MIXTURE**

The mixture is not inflammable. However, in case of fire, irritating and/or toxic fumes can be released.

**5.3 ADVICE FOR FIREFIGHTERS**

If necessary, wear self-contained breathing apparatus. Avoid that fire extinguishing water contaminates surface water and/or groundwater.

**6. ACCIDENTAL RELEASE MEASURES**

**6.1 PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES**

Keep people not involved in the intervention operations away. Ensure an adequate ventilation of the affected area. Wash with plenty of water all the body parts entered in contact with the mixture immediately. Do not handle damaged containers or the leaked product before wearing appropriate protective outfit.

**6.2 ENVIRONMENTAL PRECAUTIONS**

Avoid the contamination of surface water and soil.

**6.3 METHODS AND MATERIALS FOR CONTAINMENT AND CLEANING UP**

Absorb on inert material (sand/earth). Collect the spilled material with mechanical means and store it in suitable containers for disposal. Use water only to remove residuals, to avoid the danger of spill of product into the sewers.

**6.4 REFERENCE TO OTHER SECTIONS**

For more information on disposal see section 8 and 13.

**7. HANDLING AND STORAGE**

**7.1 PRECAUTION FOR SAFE HANDLING**

Work in well-ventilated area. Avoid the contact with eyes and skin and the inhalation of vapors.

**7.2 CONDITION FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES**

Store at 2-30°C, away from basic substances.

**7.3 SPECIFIC END USE**

None

**8. EXPOSURE CONTROL/PERSONAL PROTECTION**

**8.1 CONTROL PARAMETERS**

**Applicable exposure limits**

Description	Type	TWA/8h		STEL/15min	
		mg/m <sup>3</sup>	ppm	mg/m <sup>3</sup>	ppm
Sulphuric Acid	TLV-ACGIH 2008	0.2			

**8.2 EXPOSURE CONTROLS**

Work and handle according to the usual precautionary measures for handling chemicals. Do not eat, drink or smoke while handling the product; wash hands thoroughly with soap and water before meals and after the work shift.

**Appropriate engineering controls:**

Ensure an adequate ventilation of the working area.

**Personal Protective Equipment:**

Suggestions on the use of specific PPE are indicated below. Their choice should be made according to the use of the product and the instructions given by the supplier of the device.

Hand protections:

Chemical-resistant gloves, comply with EN 374.

Eye protections:

Safety goggles side shields

Body protections:

Work outfits for complete skin protection

Respiratory protection:

Not required under normal work activities

**8.3 ENVIRONMENTAL EXPOSURE CONTROLS**

Do not discard residuals in the environment. Discard the empty containers as hazardous waste.

**9. PHYSICAL AND CHEMICAL PROPERTIES**

**9.1 INFORMATION ON BASIC PHYSICAL AND CHEMICAL PROPERTIES**

Form	Liquid
Odour	No data available
Odour Threshold	No data available
pH	< 0.5
Freezing/Melting point	No data available
Initial boiling and melting point range	No data available
Flash point	No data available
Evaporation rate	No data available
Flammability (solid, gas)	Not applicable
Upper/lower flammability or explosive limit	Not applicable
Vapour pressure	No data available
Vapour density	No data available
Relative density	No data available
Water solubility	Soluble in water
Partition coefficient: n-octano/water	No data available
Autoignition temperature	Not applicable
Decomposition temperature	No data available
Viscosity	No data available
Explosive properties	No data available
Oxidizing properties	No data available

**9.2 OTHER SAFETY INFORMATION**

None

**10. STABILITY AND REACTIVITY**

**10.1 REACTIVITY**

No data available

**10.2 CHEMICAL STABILITY**

Stable under the recommended transport, handling and storage conditions.

**10.3 POSSIBILITY OF HAZARDOUS REACTIONS**

The mixture reacts with the basic and/or oxidizing substances with formation of heat.

**10.4 CONDITIONS TO AVOID**

Avoid contact with basic and/or oxidizing substances.

**10.5 INCOMPATIBLE MATERIALS**

Basic substances.



**MATERIAL SAFETY DATASHEET**  
**ENZY-WELL STOP SOLUTION**  
(According to Regulation (EC) No. 1907/2006)

Edition: 4  
Date: 18-Jun-2018  
Page: 6/7

**10.6 HAZARDOUS DECOMPOSITION PRODUCTS**

When heated or in case of fire, vapors that are potentially dangerous to health (SOx) may be produced.

**11. TOXICOLOGICAL INFORMATION**

**11.1 INFORMATION ON TOXICOLOGICAL EFFECTS**

The mixture is corrosive for skin and eyes.

Acute toxicity:	No data available
Skin irritation:	Corrosive for skin and mucous membranes
Eye irritation:	Severe eye damages
Subcutaneous sensitization:	No data available
Germ cell mutagenesis:	No data available
Carcinogenicity:	No data available
Reproduction toxicity:	No data available
Organ specific toxicity (on single or repeated exposure):	No data available
Acute inhalation toxicity:	No data available
Other information	No data available

**12. ECOLOGICAL INFORMATION**

**12.1 TOXICITY**

Use according to good working practices, avoiding release to the environment. In case the product reach waterways or sewers or contaminate soil or vegetation, inform the competent authorities.

**12.2 PERSISTENCE AND DEGRADABILITY**

No data available

**12.3 BIOACCUMULATIVE POTENTIAL**

No data available

**12.4 MOBILITY IN SOIL**

No data available

**12.5 RESULTS OF PBT AND VPVB ASSESSMENT**

No data available

**12.6 OTHER ADVERSE EFFECTS**

The disposal into drains or water bodies can result in a decreasing of the pH value. A low pH value harms aquatic organisms.

**13. DISPOSAL CONSIDERATIONS**

**13.1 WASTE TREATMENT METHODS**

Serum samples and reagents used should be treated as infectious waste.  
The product and its containers should be considered special waste.  
Their transport and disposal should be performed by authorized specialized companies according to the law.

**14. TRANSPORT INFORMATION**

Even if the mixture is classified as hazardous by the transport regulations (see table below), the exemptions to the classification of ADR Chapter 3.5 - *Dangerous goods packed in excepted quantities* – can be applied. These goods must be transported by vehicles authorized to the transport of hazardous goods according to the requirements of the current edition of ADR and to the national regulations. People, responsible for the loading and unloading of this dangerous good must be trained on all the risks deriving from the substance and on procedures to be taken in case of emergency.



**MATERIAL SAFETY DATASHEET**  
**ENZY-WELL STOP SOLUTION**  
(According to Regulation (EC) No. 1907/2006)

Edition: 4  
Date: 18-Jun-2018  
Page: 7/7

<b>ADR/RID:</b>			
UN N°	2796	Class:	8
Classification:	C1	Packaging group:	II
Tunnel restriction code:	(E)		
If subject to the following characteristics is ADR exempt: Exemption for limited quantities: 1L Exemption for expected quantities: <ul style="list-style-type: none"><li>Maximum quantity for the internal package: 30 ml</li><li>Maximum quantity for the external package: 500 ml</li></ul>			
<b>Sea transport (IMDG Code):</b>			
UN N°	2796	Class:	8
Packaging group:	II	EmS:	F-A, S-B
Marine polluting agent	no		
<b>Air transport (ICAO-IATA):</b>			
Class:	8	Packaging group:	II
Name:	Sulphuric Acid (not containing more than 51% of acid)		
Label:	n.8		

## 15. REGULATORY INFORMATION

### 15.1 SAFETY, HEALTH AND ENVIRONMENTAL REGULATIONS/LEGISLATION SPECIFIC FOR THE SUBSTANCE OR MIXTURE

Classification and labeling is done in accordance with the directives 1272/2008 (EC) and its subsequent amendments.

#### Legislative Decree 81/2008:

The use of this chemical entails the obligation of "Risk Assessment" by the employer in accordance with the provisions of the Decree April 9, 2008 n. 81 and subsequent amendments. Workers exposed to this chemical agent must not undergo health checks whether the results of the risk assessment show that, depending on the type and quantity of chemical agent and the method and frequency of exposure to this agent, there is only one "Low risk to the safety and irrelevant to the health" of the workers and the measures provided for in the same Decree are sufficient to reduce the risk.

### 15.2 CHEMICAL SAFETY ASSESSMENT

Not performed for the product.

## 16. OTHER INFORMATION

This product has to be used for diagnostic use only by personnel who is qualified and trained on the hazards shown in this safety sheet.

### **Text of hazard statements set out in paragraph 3**

H314 Causes severe skin burns and eye damage

#### REFERENCES:

1. Regolamento (CE) 1907/2006 del Parlamento Europeo (REACH) e s.m.i.
2. Regolamento (CE) 1272/2008 del Parlamento Europeo (CLP) e s.m.i
3. ECHA European Chemicals Agency [www.echa.europa.eu](http://www.echa.europa.eu)
4. The Merck Index
5. Handling Chemical Safety
6. Niosh - Registry of Toxic Effects of Chemical Substances
7. INRS - Fiche Toxicologique

DIESE Diagnostica Senese SpA declares that the information contained in this data sheet is based on the knowledge available to us on the date of the last version. Users must verify the suitability and accuracy of the provided information according to the specific use of the product. The document has not been accepted as a guarantee of any specific product property. Because of the use of this product is not subject to our direct control, the users, under their own responsibility, must observe the laws and regulations relating to health and safety. We accept no responsibility for improper use.

**This data sheet supersedes each previous edition.**