



QUALITY CONTROL CERTIFICATE

KIT NAME: 88000-G4 MEADOW FESCUE

Code	88000-G4		
Lot No.	852	2023-07	2022-09

REAGENTS

Label	Description	Lot	Exp	No
D0081	88000-G4 MEADOW FESCUE	852	2023-07	2X6
B0080	CALIBRATORE 88000-G4 MEADOW FESCUE	852	2023-10	1
C0080	CONTROL SERUM 88000-G4 MEADOW FESCUE	852	2023-10	1

CONTROLS


	Found results	Expected result
--	---------------	-----------------

Calibrator (mOD)	270	≥ 100
Control Serum (KU/L)	14,6	6,0-22,8

CALIBRATOR

6	7	5	8	5	2	6	5	2	7	5	7	2	2	7	7	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *SB*



CONTROL SERUM

6	7	5	8	5	2	6	3	2	6	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *SB*



The test has been carried out in accordance with the package insert

TESTED BY: Simone Balocchi *Simone Balocchi*
APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*
DATE *13.09.22*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an FDA approved method and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.