

**DIESSSE DIAGNOSTICA SENESE SpA**  
 53035 Monteriggioni (SI) – Italia – Via delle Rose 10  
 Tel +39-0577-587111 Fax +39-0577-318690



Pag 1/1  
 Sede Legale e Amministrazione  
 20144 Milano (Italy) Via A. Solari 19 sc. 6  
 Tel. +39-02-4859121 r.a. Fax +39-02-48008530

**QUALITY CONTROL CERTIFICATE**

**KIT NAME: 88000-F35 POTATO**

Code	88000-F35		
Lot No.	870	2023-09	2022-11

**REAGENTS**

Label	Description	Lot	Exp	No
D1071	88000-F35 POTATO	870	2023-09	2x6
B1070	CALIBRATORE 88000-F35 POTATO	870	2023-10	1
C1070	CONTROL SERUM 88000-F35 POTATO	870	2023-10	1


**CONTROLS**

	Found results	Expected result
Calibrator (mOD)	270	≥ 100
Control Serum (KU/L)	14,0	6,0-22,8

**CALIBRATOR**

6	6	5	8	7	0	6	5	2	7	5	7	2	2	7	7	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader **YES**  Operator Signature *SB*



**CONTROL SERUM**

6	6	5	8	7	0	6	3	2	6	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader **YES**  Operator Signature *SB*



The test has been carried out in accordance with the package insert

TESTED BY: Simone Baiocchi *Simone Baiocchi*  
 APPROVED RESPONSIBLE HEAD, Quality Control Reagents... *[Signature]*  
 DATE *22-11-2022*

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**  
 This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an FDA approved method and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.