



QUALITY CONTROL CERTIFICATE

KIT NAME 88000-F48 ONION

Code	88000-F48		
Lot No.	878	2023-10	2022-12

REAGENTS

Label	Description	Lot.	Exp.	NO
D1101	88000-F48 ONION	878	2023-10	2x6
B1100	CALIBRATORE 88000-F48 ONION	878	2024-02	1
C1100	CONTROL SERUM 88000-F48 ONION	878	2024-02	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	256	≥ 100
Control Serum (KU/L)	13,1	7,9-17,6

680878661856229300
Onion

CALIBRATOR 0.250 mL IVD

LOT 878

2024-02 2°C 8°C

DIESSA Monteriggioni Italy CE

CALIBRATOR

6	8	0	8	7	8	6	5	1	8	5	5	2	2	9	3	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *88*

CONTROL SERUM

6	8	0	8	7	8	6	2	9	5	4	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *88*

6808786629540000001
Onion

CONTROL+ 0.450 mL IVD

LOT 878

2024-02 2°C 8°C

DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *Deus*

DATE *13.12.2022*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.