

DIESE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME dsDNA-M - CHORUS

Code	86034				
Lot No.	901		2024-04		2023-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D7401	dsDNA-M	901	2024-04	6X6
B7401	CALIBRATORE dsDNA-M	901	2024-04	1
C7401	CONTROL SERUM dsDNA-M	901	2024-04	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	515	≥ 100
Control Serum (AU/mL)	58,6	38,9-86,4

CALIBRATOR

2	5	7	9	0	1	7	7	1	9	3	0	4	8	4	1	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *f*

267901771930484100  
dsDNA-M  
CALIBRATOR 0.175 mL IVD  
LOT 901  
2024-04 2°C-8°C  
DIESE Monteriggioni Italy CE

CONTROL SERUM

2	5	7	9	0	1	7	1	5	2	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *f*

267901716270000001  
dsDNA-M  
CONTROL 0.425 mL IVD  
LOT 901  
2024-04 2°C-8°C  
DIESE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Alvares*  
DATE 30/01/23

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**  
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent reagents should be handled with maximum attention.