

**QUALITY CONTROL CERTIFICATE**

**KIT NAME** HERPES SIMPLEX 1 IgG RECOMBINANT - CHORUS

Code	81023				
Lot No.	917		2024-03		2023-03

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D8681	Herpes simplex 1 IgG Recombinant	917	2024-03	6x6
B8681	CALIBRATORE Herpes simplex 1 IgG Recombinant	917	2024-03	1
C8681	CONTROL SERUM Herpes simplex 1 IgG Recombinant	917	2024-03	1

**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	434	> 100
CONTROL SERUM (INDEX)	4,0	2,7-5,0
NEGATIVE CONTROL (INDEX)	0,1	<0,7

**CALIBRATOR**

2	5	8	9	1	7	6	0	0	0	0	0	2	6	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FD*



**CONTROL SERUM**

2	5	8	9	1	7	6	2	7	5	0	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FD*



*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *FD*

DATE *15.05.23*

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**  
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.