



DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020

QUALITY CONTROL CERTIFICATE

KIT NAME MYCOPLASMA PNEUMONIAE IgM - CHORUS

Code	81035				
Lot No.	927		2024-07		2023-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D2261	MYCOPLASMA PNEUMONIAE IgM	927	2024-07	6x6
B2261	CALIBRATORE MYCOPLASMA PNEUMONIAE IgM	927	2024-07	1
C2261	CONTROL SERUM MYCOPLASMA PNEUMONIAE IgM	927	2024-07	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	510	≥ 100
Control Serum (INDEX)	1,6	1,2-2,3
Negative Control (INDEX)	0,3	≤ 0,7

CALIBRATOR

0	6	7	9	2	7	7	0	0	0	0	0	5	8	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FB*

067927700000680000  
Mycoplasma pneumoniae IgM  
CALIBRATOR 0.175 mL [IVD]  
LOT 927  
2024-07 2°C-8°C  
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	6	7	9	2	7	7	1	2	2	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FB*

067927712230000001  
Mycoplasma pneumoniae IgM  
CONTROL SERUM 0.425 mL [IVD]  
LOT 927  
2024-07 2°C-8°C  
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *FB*

DATE 26/05/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.