

DIESE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME DIESE W.W.-DAT

Code 26010

LOT.NO 154

2024-08

2023-05

REAGENTS

CODE	LABEL	DESCRIPTION	LOT	EXP	NO
V301	B1970	S. TYPHI O	158	2024-08	1
V302	E1970	S. TYPHI H	157	2024-09	1
V303	F1970	S. PARATYPHI A-O	144	2024-08	1
V304	G1970	S. PARATYPHI A-H	150	2024-08	1
V305	H1970	S. PARATYPHI B-O	145	2024-08	1
V306	I1970	S. PARATYPHI B-H	144	2024-09	1
V307	L1970	BRUCELLA	149	2024-09	1
PF21250/I	C1970	CONTROLLO POSITIVO	173	2024-08	1
PF02102/DAT	D1970	CONTROLLO NEGATIVO	137	2024-09	1

CONTROLS

	FOUND RESULTS	EXPECTED RESULTS
CONTROLLO POSITIVO	Title	Title
S. TYPHI O	1/320	1/320 +/- 1 Dil
S. TYPHI H	1/160	1/160 +/- 1 Dil
S. PARATYPHI A-O	1/320	1/320 +/- 1 Dil
S. PARATYPHI A-H	1/320	1/320 +/- 1 Dil
S. PARATYPHI B-O	1/320	1/320 +/- 1 Dil
S. PARATYPHI B-H	1/320	1/320 +/- 1 Dil
BRUCELLA	1/320	1/320 +/- 1 Dil

*Evaluation of the method according to the instruction for use.  
This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED BY: RESPONSIBLE HEAD, Quality Control Reagents

DATE: 07/07/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.