

QUALITY CONTROL CERTIFICATE

KIT NAME 88000-D2 DERMATOPHAGOIDES FARINAE

Code	88000-D2		
Lot No.	943	2024-06	2023-08

REAGENTS

Label	Description	Lot.	Exp.	NO
D9801	88000-D2 DERMATOPHAGOIDES FARINAE	943	2024-06	2X6
B9800	CALIBRATORE 88000-D2 DERMATOPHAGOIDES FARINAE	943	2024-06	1
C9800	CONTROL SERUM 88000-D2 DERMATOPHAGOIDES FARINAE	943	2024-06	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	290	≥ 100
Control Serum (KU/L)	17,2	10,0-41,5

CALIBRATOR

6 0 3 9 4 3 6 5 2 1 6 6 1 7 8 1 0 0

Label checked with barcode reader YES Operator Signature *SB*



CONTROL SERUM

6 0 3 9 4 3 6 4 1 7 7 0 0 0 0 0 0 1

Label checked with barcode reader YES Operator Signature *SB*



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *SB*

DATE *01.08.2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.