



QUALITY CONTROL CERTIFICATE

KIT NAME a-TPO - CHORUS

Code	86074/12		
Lot No.	946	2024-06	2023-06

REAGENTS

Label	Description	Lot.	Exp.	NO
D3591	a-TPO	946	2024-06	2x6
B3591	CALIBRATORE a-TPO	946	2024-06	1
C3591	CONTROL SERUM a-TPO	946	2024-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	379	≥ 100
Control Serum (IU/mL)	310,3	233,3-533,3

CALIBRATOR

3	0	1	9	4	6	6	5	1	9	3	7	3	4	5	2	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES  Operator Signature *SB*



CONTROL SERUM

3	0	1	9	4	6	6	2	1	3	8	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *SB*



**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *SB*

DATE *14.09.2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.