

QUALITY CONTROL CERTIFICATE

KIT NAME ANA-Screen - CHORUS

Code	86014		
Lot No.	947	2024-06	2023-06

REAGENTS

Label	Description	Lot.	Exp.	NO
D3141	ANA-Screen	947	2024-06	6X6
B3141	CALIBRATORE ANA-Screen	947	2024-06	1
C3141	CONTROL SERUM ANA-Screen	947	2024-06	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	330	≥ 100
Control Serum (Ndx)	3,5	2,3-4,2

CALIBRATOR

1	1	5	9	4	7	6	0	0	0	0	0	0	3	1	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature SB

11594760000310000  
ANA-Screen  
CALIBRATOR 0.175 mL IVD  
LOT 947  
2024-06 2°C-8°C  
DIESSE Monteriggioni Italy CE

CONTROL SERUM

1	1	5	9	4	7	6	2	3	4	2	0	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature SB

115947623420000001  
ANA-Screen  
CONTROL+ 0.425 mL IVD  
LOT 947  
2024-06 2°C-8°C  
DIESSE Monteriggioni Italy CE

**This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents SB

DATE 15.09.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.