



DIESSSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020

QUALITY CONTROL CERTIFICATE

KIT NAME 88000-F76 ALPHA-LACTALBULMIN

Code	88000-F76		
Lot No.	963	2024-08	2023-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D9861	88000-F76 ALPHA-LACTALBULMIN	963	2024-08	2X6
B9860	CALIBRATORE 88000-F76 ALPHA-LACTALBULMIN	963	2024-10	1
C9860	CONTROL SERUM 88000-F76 ALPHA-LACTALBULMIN	963	2024-10	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	242	≥ 100
Control Serum (KU/L)	12,1	7,5-15,7

CALIBRATOR

6	1	9	9	6	3	6	5	1	8	6	2	2	4	8	7	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FB*

619963651862248700  
Alpha-lactalbumin  
CALIBRATOR 0.250 mL [I]VD  
LOT 963  
2024-10 2°C-8°C  
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

6	1	9	9	6	3	6	3	0	5	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FB*

619963630550000001  
Alpha-lactalbumin  
CONTROL 0.450 mL [I]VD  
LOT 963  
2024-10 2°C-8°C  
DIESSSE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Bene*

DATE 03/10/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.