



QUALITY CONTROL CERTIFICATE

KIT NAME VARICELLA IgM CHORUS

Code	81079			
Lot No.	955	2024-10	2023-07	

REAGENTS

Label	Description	Lot.	Exp.	NO
D2481	VARICELLA IgM	955	2024-10	6x6
B2481	CALIBRATORE VARICELLA IgM	955	2024-10	1
C2481	CONTROL SERUM VARICELLA IgM	955	2024-10	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	615	≥ 100
Control Serum (Index)	2,4	1,7-3,1
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

1	3	4	9	5	5	7	0	0	0	0	0	0	4	2	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

134955700000420000
Varicella IgM
CALIBRATOR 0.175 mL [IVD]
LOT 955
2024-10 2°C-8°C
DIESE Monteriggioni Italy CE

CONTROL SERUM

1	3	4	9	5	5	7	1	7	3	1	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

134955717310000001
Varicella IgM
CONTROL+ 0.425 mL [IVD]
LOT 955
2024-10 2°C-8°C
DIESE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents Bene
DATE 12/10/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive