



QUALITY CONTROL CERTIFICATE

KIT NAME TOXOPLASMA IgG - CHORUS

Code	81040		
Lot No.	971	2024-09	2023-09

REAGENTS

Label	Description	Lot.	Exp.	NO
D2271	TOXOPLASMA IgG	971	2024-09	6X6
B2272	CALIBRATORE TOXOPLASMA IgG	971	2024-09	1
C2271	CONTROL SERUM TOXOPLASMA IgG	971	2024-09	1

CONTROLS


	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	267	≥ 100
Control Serum (IU/mL)	35,2	20,0-40,0
Negative Control (IU/mL)	<5	<5

CALIBRATOR

0	8	2	9	7	1	6	3	2	5	4	5	3	5	7	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *HP*



CONTROL SERUM

0	8	2	9	7	1	6	2	0	3	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *[Signature]*
DATE *13/10/23*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive