

QUALITY CONTROL CERTIFICATE

KIT NAME HELICOBACTER PYLORI IgG -CHORUS

Code	81060
Lot No.	978
	2024-10
	2023-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D2391	Helicobacter Pylori IgG	978	2024-10	6X6
B2392	CALIBRATORE Helicobacter Pylori IgG	978	2024-10	1
C2391	CONTROL SERUM Helicobacter Pylori IgG	978	2024-10	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	310	≥ 100
Control Serum (UA/mL)	78,0	55,3-141,7
Negative Control (UA/mL)	<5,0	< 5

CALIBRATOR

0	0	8	9	7	8	6	1	2	3	4	2	2	8	5	4	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *Mb*

008978612342285400
Helicobacter pylori IgG
CALIBRATOR 0.425 mL IVD
LOT 978
2024-10 2°C-8°C
DIESSA Monteriggioni Italy CE

CONTROL SERUM

0	0	8	9	7	8	6	2	5	4	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *Mb*

008978625470000001
Helicobacter pylori IgG
CONTROL+ 0.425 mL IVD
LOT 978
2024-10 2°C-8°C
DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Mb*

DATE 23.10.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.