



QUALITY CONTROL CERTIFICATE

KIT NAME ADENOVIRUS IgA - CHORUS

Code	81198		
Lot No.	978	2024-10	2023-10

REAGENTS

Label	Description	Lot.	Exp.	NO
D7611	Adenovirus IgA	978	2024-10	6X6
B7611	CALIBRATORE Adenovirus IgA	978	2024-10	1
C7611	CONTROL SERUM Adenovirus IgA	978	2024-10	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	435	≥ 100
Control Serum (Ndx)	2,6	1,8-3,3
Negative Control (Ndx)	0,1	<0,7

CALIBRATOR

0	7	5	9	7	8	6	0	0	0	0	0	4	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

07597860000400000
Adenovirus IgA
CALIBRATOR 0.175 mL IVD
LOT 978
 2024-10 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	7	5	9	7	8	6	1	8	3	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

07597861833000001
Adenovirus IgA
CONTROL 0.425 mL IVD
LOT 978
 2024-10 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.

DATE 07-11-2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.