



DIESSSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME CARDIOLIPIN-G - CHORUS

Code	86046		
Lot No.	987	2024-11	2023-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D3311	CARDIOLIPIN-G	987	2024-11	6X6
B3311	CALIBRATORE CARDIOLIPIN-G	987	2024-11	1
C3311	CONTROL SERUM CARDIOLIPIN-G	987	2024-11	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	351	≥ 100
Control Serum (GPLU/mL)	46,0	34,3-85,3

CALIBRATOR

3	3	5	9	8	7	6	4	2	6	4	0	2	5	5	9	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature RB

335987642640255900
Cardiolipin-G
CALIBRATOR 0.175 mL IVD
LOT 987
2024-11 2°C-8°C
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

3	3	5	9	8	7	6	2	8	5	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature RB

335987628520000001
Cardiolipin-G
CONTROL+ 0.425 mL IVD
LOT 987
2024-11 2°C-8°C
DIESSSE Monteriggioni Italy CE

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. lesca

DATE 22/11/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.