

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME HERPES SIMPLEX 1+2 IgM -CHORUS

Code	81021		
Lot No.	987	2024-11	2023-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D2221	Herpes simplex 1+2 IgM	987	2024-11	6x6
B2221	CALIBRATORE Herpes simplex 1+2 IgM	987	2024-11	1
C2221	CONTROL SERUM Herpes simplex 1+2 IgM	987	2024-11	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	590	≥ 100
Control Serum (INDEX)	2,9	2,0-3,6
Negative Control (INDEX)	0,0	<0,7

CALIBRATOR

1	5	6	9	8	7	6	0	0	0	0	0	0	3	6	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature HB

156987600000360000
Herpes simplex 1+2 IgM
CALIBRATOR 0.175 mL [IVD]
LOT 987
2024-11 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

1	5	6	9	8	7	6	2	0	3	6	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature HB

15698762036000001
Herpes simplex 1+2 IgM
CONTROL+ 0.425 mL [IVD]
LOT 987
2024-11 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... HB

DATE 30.11.2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.