



DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME Epstein-Barr VCA IgG - CHORUS

Code	81055		
Lot No.	987	2025-02	2023-11

REAGENTS

Label	Description	Lot.	Exp.	NO
D2341	Epstein-Barr VCA IgG	987	2025-02	6X6
B2342	CALIBRATORE Epstein-Barr VCA IgG	987	2025-02	1
C2341	CONTROL SERUM Epstein-Barr VCA IgG	987	2025-02	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	350	≥ 100
Control Serum (INDEX)	4,5	3,3-6,2
Negative Control (INDEX)	0,0	<0,7

CALIBRATOR

0	0	5	9	8	7	7	0	0	0	0	0	2	1	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

005987700000210000
Epstein-Barr VCA IgG
CALIBRATOR 0.425 mL IVD
LOT 987
2025-02 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	0	5	9	8	7	7	3	3	6	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature MB

005987733620000001
Epstein-Barr VCA IgG
CONTROL+ 0.425 mL IVD
LOT 987
2025-02 2°C-8°C
DIESSE Monteriggioni Italy CE

**This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents. [Signature]

DATE 30/11/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.