



QUALITY CONTROL CERTIFICATE

KIT NAME **CYTOMEGALOVIRUS IgG -CHORUS**

Code	81010		
Lot No.	996	2024-12	2023-12

REAGENTS

Label	Description	Lot.	Exp.	NO
D2161	Cytomegalovirus IgG	996	2024-12	6x6
B2162	CALIBRATORE Cytomegalovirus IgG	996	2024-12	1
C2161	CONTROL SERUM Cytomegalovirus IgG	996	2024-12	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	310	≥ 100
Control Serum (IU/mL)	15,5	7,7-20,0
Negative Control (IU/mL)	<0,50	<0,5

CALIBRATOR

0	0	1	9	9	6	6	2	2	4	3	5	2	4	5	5	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature FB

001996622435245500  
Cytomegalovirus IgG  
CALIBRATOR 0.425 mL IVD  
LOT 996  
 2024-12 2°C-8°C  
DIESSE Monteriggioni Italy CE 0123 B2162

CONTROL SERUM

0	0	1	9	9	6	6	3	0	5	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature FB

001996630550000001  
Cytomegalovirus IgG  
CONTROL+ 0.425 mL IVD  
LOT 996  
 2024-12 2°C-8°C  
DIESSE Monteriggioni Italy CE 0123 C2161

*This batch has been manufactured and checked according to quality requirements.  
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents FB

DATE 11-01-2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.