

DIESSE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale  
20144 Milano (Italy) Via A. Solari 19 sc. 6

**QUALITY CONTROL CERTIFICATE**

**KIT NAME** HERPES SIMPLEX1 IgG RECOMBINANT - CHORUS

Code	81023		
Lot No.	996	2024-12	2023-12

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D8681	Herpes simplex 1 IgG Recombinant	996	2024-12	6x6
B8681	CALIBRATORE Herpes simplex 1 IgG Recombinant	996	2024-12	1
C8681	CONTROL SERUM Herpes simplex 1 IgG Recombinant	996	2024-12	1


**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	456	> 100
CONTROL SERUM (INDEX)	2,4	1,7-3,1
NEGATIVE CONTROL (INDEX)	0,1	<0,7

**CALIBRATOR**

2	5	8	9	9	6	6	0	0	0	0	0	4	3	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES  Operator Signature FB



**CONTROL SERUM**

2	5	8	9	9	6	6	1	7	3	1	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature FB



**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 05/01/2024

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**  
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.