



QUALITY CONTROL CERTIFICATE

KIT NAME HERPES SIMPLEX 2 IgG RECOMBINANT - CHORUS

Code	81029		
Lot No.	001	2025-01	2024-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D8690	Herpes simplex 2 IgG Recombinant	001	2025-01	6X6
B8690	CALIBRATORE Herpes simplex 2 IgG Recombinant	001	2025-01	1
C8690	CONTROL SERUM Herpes simplex 2 IgG Recombinant	001	2025-01	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	490	≥ 100
Control Serum (INDEX)	1,8	1,4-2,7
Negative Control (INDEX)	0,0	<0,7

CALIBRATOR

2	6	0	0	0	1	6	0	0	0	0	0	4	9	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

26000160000490000  
Herpes simplex 2 IgG Recombinant  
CALIBRATOR 0.175 mL IVD  
LOT 001  
2025-01 2°C-8°C  
DIESSA Monteriggioni Italy CE

CONTROL SERUM

2	6	0	0	0	1	6	1	4	2	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature MB

260001614270000001  
Herpes simplex 2 IgG Recombinant  
CONTROL+ 0.425 mL IVD  
LOT 001  
2025-01 2°C-8°C  
DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents MB

DATE 14/02/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.