

**QUALITY CONTROL CERTIFICATE**

**KIT NAME** 88000-W8 DANDELION

Code	88000-W8		
Lot No.	002	2025-01	2024-02

**REAGENTS**

Label	Description	Lot.	Exp.	NO
D6900	88000-W8 DANDELION	002	2025-01	2X6
B6900	CALIBRATORE 88000-W8 DANDELION	002	2025-02	1
C6900	CONTROL SERUM 88000-W8 DANDELION	002	2025-02	1


**CONTROLS**

	Found results	Expected results
Calibrator (mOD)	284	≥ 100
Control Serum (KU/L)	8,2	5,8-14,3

**CALIBRATOR**

7	9	9	0	0	2	6	5	2	1	6	0	2	6	8	1	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---


Label checked with barcode reader YES  Operator Signature ..... *SS* .....



**CONTROL SERUM**

7	9	9	0	0	2	6	2	7	5	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature ..... *SS* .....



**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *SS* .....

DATE *15.02.2024*

**WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!**

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.