

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME ENZYWELL RUBELLA IgG AVIDITY

Code	91095		
LOT.NO	080-A	2025-04	2024-01

REAGENTS

CODE	LABEL	DESCRIPTION	LOT	EXP	NO
PF93030	007580	MICROPLATE strips	205	2025-10	3
V140	007570	CONJUGATE	245	2025-04	1
PF91094	007550	CONTROL	138	2025-07	1
PF91091	007590	AVIDITY BUFFER	150	2025-07	1
PF93619	1034	SUBSTRATE	383	2025-05	1
PF93611	009660	DILUENT 2	232	2025-05	1
PF93603	1035	WASHING BUFFER 10x	468	2025-06	1
PF93602	1032/RISK0	STOP SOLUTION	218-6	2025-06	1

CONTROLS

	FOUND RESULTS	EXPECTED RESULTS
CONTROL	56%	> 45%

*Evaluatin of the method according to the package inset in the attached document
This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED BY: RESPONSIBLE HEAD, Quality Control Reagents

DATE: 15.02.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.



CERTIFICATE OF ANALYSIS

KIT:	RUBELLA IgG AVIDITY	Code	91095
LOT:	080		
DATE OF TEST:	17/01/2024		

METHOD: ACCORDING TO THE PACKAGE INSERT

	<i>FOUND</i>	<i>EXPECTED</i>
--	---------------------	------------------------

1. CONTROL OF COMPLETE KIT


High Avidity Control	56%	> 45%
----------------------	------------	-------

Good function with high and low avidity sera and Internal Standards, in accordance with specifications

<input checked="" type="checkbox"/>	SATISFACTORY	<input type="checkbox"/>	UNSATISFACTORY
-------------------------------------	--------------	--------------------------	----------------

FINAL RESULT: ACCEPTED

Date: 17. d. 2024

Signed: 
Released: Q.C.Manager