

DIESE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME DIESE W.W.-DAT

Code	26010		
LOT.NO	155	2025-05	2024-02

REAGENTS

CODE	LABEL	DESCRIPTION	LOT	EXP	NO
V301	B1970	S. TYPHI O	159	2025-05	1
V302	E1970	S. TYPHI H	158	2025-05	1
V303	F1970	S. PARATYPHI A-O	145	2025-05	1
V304	G1970	S. PARATYPHI A-H	151	2025-05	1
V305	H1970	S. PARATYPHI B-O	146	2025-05	1
V306	I1970	S. PARATYPHI B-H	145	2025-05	1
V307	L1970	BRUCELLA	150	2025-05	1
PF21250/1	C1970	CONTROLLO POSITIVO	174	2025-05	1
PF02102/DAT	D1970	CONTROLLO NEGATIVO	138	2025-05	1

CONTROLS

	FOUND RESULTS	EXPECTED RESULTS
CONTROLLO POSITIVO	Title	Title
S. TYPHI O	1/160	1/160 +/- 1 Dil
S. TYPHI H	1/320	1/320 +/- 1 Dil
S. PARATYPHI A-O	1/320	1/320 +/- 1 Dil
S. PARATYPHI A-H	1/320	1/320 +/- 1 Dil
S. PARATYPHI B-O	1/320	1/320 +/- 1 Dil
S. PARATYPHI B-H	1/320	1/320 +/- 1 Dil
BRUCELLA	1/320	1/320 +/- 1 Dil

*Evaluation of the method according to the instruction for use.
This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED BY: RESPONSIBLE HEAD, Quality Control Reagents

DATE: 13/02/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.