



QUALITY CONTROL CERTIFICATE

KIT NAME **Cenp-B - CHORUS**

Code	86028		
Lot No.	007	2025-02	2024-02

REAGENTS

Label	Description	Lot.	Exp.	NO
D3241	Cenp-B	007	2025-02	6X6
B3241	CALIBRATORE Cenp-B	007	2025-02	1
C3241	CONTROL SERUM Cenp-B	007	2025-02	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	420	≥ 100
Control Serum (AU/mL)	65,5	33,3-96,4

CALIBRATOR

1	9	6	0	0	7	6	1	2	0	2	9	3	3	4	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*

196007612029334000
Cenp-B
CALIBRATOR 0.176 mL IVD
LOT 007
 2025-02 B0241
DIESSE Monteriggioni Italy

CONTROL SERUM

1	9	6	0	0	7	6	2	1	3	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*

196007621390000001
Cenp-B
CONTROL+ 0.425 mL IVD
LOT 007
 2025-02 C0241
DIESSE Monteriggioni Italy

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE *16.02.2024*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.