



DIESSSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

QUALITY CONTROL CERTIFICATE

KIT NAME RUBELLA IgM - CHORUS

Code	81031		
Lot No.	002	2025-01	2024-01

REAGENTS

Label	Description	Lot.	Exp.	NO
D2241	RUBELLA IgM	002	2025-01	6X6
B2241	CALIBRATORE RUBELLA IgM	002	2025-01	1
C2241	CONTROL SERUM	002	2025-01	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	590	≥ 100
Control Serum (INDEX)	2,3	1,7-3,2
Negative Control (INDEX)	0,0	< 0,7

CALIBRATOR

1	3	5	0	0	2	6	0	0	0	0	0	4	1	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature SB

13500260000410000
Rubella IgM
CALIBRATOR 0.175 mL [IVD]
LOT 002
2025-01 2°C 8°C
DIESSSE Monteriggioni Italy CE 0123 B2241

CONTROL SERUM

1	3	5	0	0	2	6	1	7	3	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature SB

135002617320000001
Rubella IgM
CONTROL+ 0.425 mL [IVD]
LOT 002
2025-01 2°C 8°C
DIESSSE Monteriggioni Italy CE 0123 C2241

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents SB

DATE 23.02.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.