



DIESSSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code	86600			
Lot No.	010	2025-02	2024-02	

REAGENTS

Label	Description	Lot.	Exp.	NO
D9500	CALPROTECTIN	010	2025-02	2 X 6
B9500	CALIBRATORE CALPROTECTIN	010	2025-02	1
C9500	CONTROL CALPROTECTIN	010	2025-02	1
G9502	SAMPLE DILUENT	010	2026-01	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	390	≥ 100
Control (mg/kg)	164,7	100,0-225,0

CALIBRATOR

9	0	0	0	1	0	6	2	1	7	2	3	4	1	3	5	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *FB*

900010621723413500
Calprotectin
CALIBRATOR 0.300 mL IVD
LOT 010
2025-02
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

9	0	0	0	1	0	6	1	7	3	2	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *FB*

900010617320000001
Calprotectin
CONTROL+ 0.600 mL IVD
LOT 010
2025-02
DIESSSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *[Signature]*

DATE *10/04/2024*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.