



QUALITY CONTROL CERTIFICATE

KIT NAME Scl-70 - CHORUS

Code	86026		
Lot No.	025	2025-04	2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D3221	Scl-70	025	2025-04	6X6
B3221	CALIBRATORE Scl-70	025	2025-04	1
C3221	CONTROL SERUM Scl-70	025	2025-04	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	378	≥ 100
Control Serum (AU/mL)	37,2	27,0-70,9

CALIBRATOR

1	9	5	0	2	5	6	3	2	0	3	0	4	2	4	1	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

195025632030424100
Scl-70
CALIBRATOR 0.175 mL [IVD]
LOT 025
 2025-04 2°C-8°C
DIESSE Monteriggioni Italy CE
B3221

CONTROL SERUM

1	9	5	0	2	5	6	1	8	3	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

195025618330000001
Scl-70
CONTROL+ 0.425 mL [IVD]
LOT 025
 2025-04 2°C-8°C
DIESSE Monteriggioni Italy CE
C3221

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 13/05/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.