

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

DIESSE

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code 86600

Lot No. 011

2025-02

2024-02

REAGENTS

Label	Description	Lot.	Exp.	NO
D9500	CALPROTECTIN	011	2025-02	2 X 6
B9500	CALIBRATORE CALPROTECTIN	011	2025-02	1
C9500	CONTROL CALPROTECTIN	011	2025-02	1
G9502	SAMPLE DILUENT	011	2026-02	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	474	≥ 100
Control (mg/kg)	130,3	81,2-195,5

CALIBRATOR

9 0 0 0 1 1 6 4 1 8 2 3 4 7 3 4 0 0

Label checked with barcode reader YES ☒ Operator Signature MIB



CONTROL SERUM

9 0 0 0 1 1 6 1 5 2 8 0 0 0 0 0 0 1

Label checked with barcode reader YES ☒ Operator Signature HB



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 10.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.