

QUALITY CONTROL CERTIFICATE

KIT NAME CYTOMEGALOVIRUS IgM CAPTURE - CHORUS

Code 81012

Lot No. 024



2025-04



2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D3920	CYTOMEGALOVIRUS IgM CAPTURE	024	2025-04	6X6
B3920	CALIBRATORE CYTOMEGALOVIRUS IgM CAPTURE	024	2025-04	1
C3920	CONTROL SERUM CYTOMEGALOVIRUS IgM CAPTURE	024	2025-04	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	400	≥ 100
Control Serum (Index)	2,7	1,9-3,6

CALIBRATOR

9 7 1 0 2 4 6 0 0 0 0 0 3 7 0 0 0 0

Label checked with barcode reader

YES



Operator Signature FB



CONTROL SERUM

9 7 1 0 2 4 6 1 9 3 6 0 0 0 0 0 0 1

Label checked with barcode reader

YES



Operator Signature FB



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents... FB

DATE 13/05/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.