

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME **HELICOBACTER PYLORI IgA-CHORUS**

Code	81062		
Lot No.	026	2025-04	2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D2402	Helicobacter Pylori IgA	026	2025-04	6X6
B2402	CALIBRATORE Helicobacter Pylori IgA	026	2025-04	1
C2402	CONTROL SERUM Helicobacter Pylori IgA	026	2025-04	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	250	≥ 100
Control Serum (INDEX)	4,1	2,4-4,5
Negative Control (INDEX)	0,0	$< 0,7$

CALIBRATOR

0	2	7	0	2	6	6	0	0	0	0	0	0	2	9	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature

02702660000290000
HELICOBACTER PYLORI IgA
CALIBRATOR 0.175 mL IVD
LOT 026
 2025-04 2°C 8°C
DIESSE Monteriggioni Italy CE 0123
B2402

CONTROL SERUM

0	2	7	0	2	6	6	2	4	4	5	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES ☒ Operator Signature

02702662445000001
HELICOBACTER PYLORI IgA
CONTROL 0.425 mL IVD
LOT 026
 2025-04 2°C 4°C
DIESSE Monteriggioni Italy CE 0123
C2402

*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 20.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.