

QUALITY CONTROL CERTIFICATE

KIT NAME TOXOPLASMA IgM - CHORUS

Code	81041		
Lot No.	017	2025-03	2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D2281	TOXOPLASMA IgM	017	2025-03	6X6
B2281	CALIBRATORE TOXOPLASMA IgM	017	2025-03	1
C2281	CONTROL SERUM TOXOPLASMA IgM	017	2025-03	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	465	≥ 100
Control Serum (INDEX)	1,9	1,5-2,8
Negative Control (INDEX)	0,1	<0,7

CALIBRATOR

1	3	1	0	1	7	6	0	0	0	0	0	0	4	6	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

131017600000460000
Toxoplasma IgM
CALIBRATOR 0.175 mL IVD
LOT 017
2025-03 2°C 8°C
DIESSE Monteriggioni Italy CE 0123 B2281

CONTROL SERUM

1	3	1	0	1	7	6	1	5	2	8	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

131017615280000001
Toxoplasma IgM
CONTROL 0.425 mL IVD
LOT 017
2025-03 2°C 8°C
DIESSE Monteriggioni Italy CE 0123 C2281

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 21.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.