

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020



Sede Legale e Amministrazione
20144 Milano (Italy) Via A. Solari 19 sc. 6

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME ENZYWELL RUBELLA IgG

Code 91030

LOT.NO 286



2025-07



2024-05

REAGENTS

CODE	LABEL	DESCRIPTION	LOT	EXP	NO
PF93030	006510	MICROPLATE	207	2026-04	1
PF93507	006500	CONJUGATE	246	2025-08	1
PF93910	1036	IgG NEGATIVE CONTROL	180	2025-11	1
PF93831	006520	CALIBRATOR 1	203	2025-11	1
PF91830	006540	CALIBRATOR 2	203	2025-11	1
PF93833	006550	CALIBRATOR 3	203	2025-11	1
PF93834	006560	CALIBRATOR 4	203	2025-11	1
PF93835	006530	CALIBRATOR 5	203	2025-11	1
PF93619	1034	SUBSTRATE-HS	385	2025-08	1
PF93611	009660	DILUENT 2	233	2025-07	1
PF93603	1035	WASHING BUFFER 10x	473	2025-08	1
PF93602	1032/RISKO	STOP SOLUTION	220	2026-03	1

CONTROLS

IU/ml	FOUND RESULTS	EXPECTED RESULTS
IgG NEG CONTROL 0	0.003	≤ 0,6 CUT-OFF
CALIBRATOR 1 5	0.268	< CUT-OFF
CALIBRATOR 2 10 (Cut Off)	0.555	≥ 0.200
CALIBRATOR 3 50	1.126	
CALIBRATOR 4 100	1.710	
CALIBRATOR 5 200	2.456	≥ 1.500
CV %	12%	≤ 15%

*Evaluation of the method according to the package inset in the attached document
This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED BY: RESPONSIBLE HEAD, Quality Control Reagents

DATE: 21.05.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested by an and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.

IO 09 - 010-V



CODE 91030

CODE 91176

QUALITY CONTROL CERTIFICATE

Annex: Method Evaluation

NAME KIT RUBELLA IgG
LOT. 286
DATE OF TEST 17/05/2024

METHOD: ACCORDING TO THE PACKAGE INSERT

1. SOLID PHASE	FOUND	EXPECTED
A) UNIFORMITY BETWEEN PLATES		
1st PLATE	OD450 nm 2.483 CV% 5%	< 20%
2nd PLATE	OD 450 nm 2.613 CV% 6%	< 20%
B) UNIFORMITY WITHIN PLATE	OD 450 nm 2.548 CV% 6%	< 15%
2. BULK CONTROL date: 14/05/2024		
A) Average OD of Cut Off Calibrator	OD 450 nm 0.582	≥ 0.200
B) Index of Negative and Positive Sera	accepted	According to specification
C) Sensitivity with working calibrators	accepted	According to specification
3. CONTROL OF COMPLETE KIT date: 17/05/2024		
A) Substrate Blank	OD 450 nm 0.046	≤ 0.100
Calibrator 0 IU/ml	OD 450 nm 0.003	≤ 0,6 Cut-off
Calibrator 5 IU/ml	OD 450 nm 0.268	< Cut-off
Calibrator 10 IU/ml (Cut-off)	OD 450 nm 0.555	≥ 0.200
Calibrator 50 IU/ml	OD 450 nm 1.126	
Calibrator 100 IU/ml	OD 450 nm 1.710	
Calibrator 200 IU/ml	OD 450 nm 2.456	≥ 1.500
B) HOMOGENEITY OF CUT-OFF	No. 6 OD 450 nm 0.602 CV% 12%	≤ 15%
C) PRECISION OF FILLING PROCESS		



SATISFACTORY



UNSATISFACTORY

FINAL RESULT: ACCEPTED

Date: 20.05.2024

Signed:

Released: Q.C. Manager