

DIESSSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME MPO - CHORUS

Code	86068		
Lot No.	033	2025-05	2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D3532	MPO	033	2025-05	6X6
B3532	CALIBRATORE MPO	033	2025-05	1
C3532	CONTROL SERUM MPO	033	2025-05	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	570	≥ 100
Control Serum (AU/mL)	35,7	21,7-53,3

CALIBRATOR

2	5	3	0	3	3	6	1	1	9	2	8	5	1	4	1	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*

253033611926514100
MPO
CALIBRATOR 0,176 mL [IVD]
LOT 033
2025-05 2°C-8°C
DIESSSE Monteriggioni Italy CE 0123

CONTROL SERUM

2	5	3	0	3	3	6	1	4	2	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *MB*

253033614260000001
MPO
CONTROL [+] 0,425 mL [IVD]
LOT 033
2025-05 2°C-8°C
DIESSSE Monteriggioni Italy CE 0123

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*

DATE *28.05.2024*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.