



QUALITY CONTROL CERTIFICATE

KIT NAME 88000-F48 ONION

Code	88000-F48		
Lot No.	036	2025-05	2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D1101	88000-F48 ONION	036	2025-05	2x6
B1100	CALIBRATORE 88000-F48 ONION	036	2025-06	1
C1100	CONTROL SERUM 88000-F48 ONION	036	2025-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	300	≥ 100
Control Serum (KU/L)	10,7	6,6-15,0

CALIBRATOR

6 8 0 0 3 6 6 5 1 9 5 8 2 7 7 8 0 0

Label checked with barcode reader YES Operator Signature *FB*

680036651958277800
Onion
CALIBRATOR 0.280 mL IVD
LOT 036
2025-06 2°C-8°C
DIESSA Monteriggioni Italy CE

CONTROL SERUM

6 8 0 0 3 6 6 2 7 5 0 0 0 0 0 0 0 1

Label checked with barcode reader YES Operator Signature *FB*

680036627500000001
Onion
CONTROL 0.450 mL IVD
LOT 036
2025-06 2°C-8°C
DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents *[Signature]*
DATE 31/05/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.