



QUALITY CONTROL CERTIFICATE

KIT NAME TOSCANA VIRUS IgM - CHORUS

Code	81068		
Lot No.	033	2025-05	2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D2422	TOSCANA VIRUS IgM	033	2025-05	6X6
B2422	CALIBRATORE TOSCANA VIRUS IgM	033	2025-05	2
C2422	CONTROL SERUM TOSCANA VIRUS IgM	033	2025-05	2

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	428	≥ 100
Control Serum (Index)	3,1	1,9-3,6
Negative Control (Index)	0,1	< 0,7

CALIBRATOR

0	2	5	0	3	3	6	0	0	0	0	0	3	7	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature FB



CONTROL SERUM

0	2	5	0	3	3	6	1	9	3	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature FB



This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents FB

DATE 10/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.