

DIESSA DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME CHORUS Aspergillus Galactomannan Ag

Code	81178		
Lot No.	033	2025-05	2024-06

REAGENTS

Label	Description	Lot	Exp.	NO
D4070	ASPERGILLUS Ag	033	2025-05	6X6
B4070	CALIBRATORE ASPERGILLUS Ag	033	2025-05	1
C4070	CONTROL SERUM ASPERGILLUS Ag	033	2025-05	2
G4070	SAMPLE DILUENT ASPERGILLUS Ag	LP2402	2025-09	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	770	≥ 100
Control Serum (Index)	1,8	1,2-2,3

CALIBRATOR

9	7	7	0	3	3	6	0	0	0	0	0	5	7	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature FB

97703360000570000
Aspergillus
Galactomannan Ag
CALIBRATOR 0.950 mL [IVD]
LOT 033
2025-05
DIESSA Monteriggioni Italy CE

CONTROL SERUM

9	7	7	0	3	3	6	1	2	2	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature FB

977033612230000001
Aspergillus
Galactomannan Ag
CONTROL 1.250 mL [IVD]
LOT 033
2025-05
DIESSA Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents [Signature]

DATE 14/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.