

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554

DIESSE

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME CALPROTECTIN - CHORUS

Code 86600

Lot No. 034

2025-05

2024-05

REAGENTS

Label	Description	Lot.	Exp.	NO
D9500	CALPROTECTIN	034	2025-05	2 X 6
B9500	CALIBRATORE CALPROTECTIN	034	2025-10	1
C9500	CONTROL CALPROTECTIN	034	2025-10	1
G9502	SAMPLE DILUENT	034	2026-05	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	356	≥ 100
Control (mg/kg)	170,8	100,0-250,0

CALIBRATOR

9	0	0	0	3	4	6	2	1	8	2	3	3	9	3	3	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES

Operator Signature




CONTROL SERUM

9	0	0	0	3	4	6	1	8	3	3	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader

YES

Operator Signature



This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 18.06.2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.