

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME **HELICOBACTER PYLORI IgA-CHORUS**

Code **81062**

Lot No. **028**



2025-04



2024-04

REAGENTS

Label	Description	Lot.	Exp.	NO
D2402	Helicobacter Pylori IgA	028	2025-04	6X6
B2402	CALIBRATORE Helicobacter Pylori IgA	028	2025-04	1
C2402	CONTROL SERUM Helicobacter Pylori IgA	028	2025-04	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	310	≥ 100
Control Serum (INDEX)	3,0	2,1-3,9
Negative Control (INDEX)	0,0	$< 0,7$

CALIBRATOR

0 2 7 0 2 8 6 0 0 0 0 0 3 4 0 0 0 0

Label checked with barcode reader

YES



Operator Signature

MB



CONTROL SERUM

0 2 7 0 2 8 6 2 1 3 9 0 0 0 0 0 0 1

Label checked with barcode reader

YES



Operator Signature

MB



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

dfme

DATE

20/06/2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.