

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554



QUALITY CONTROL CERTIFICATE

KIT NAME VARICELLA IgM CHORUS

Code	81079		
Lot No.	034	2025-08	2024-05

REAGENTS 034

Label	Description	Lot	Exp.	NO
D2481	VARICELLA IgM	034	2025-08	6x6
B2481	CALIBRATORE VARICELLA IgM	034	2025-08	1
C2481	CONTROL SERUM VARICELLA IgM	034	2025-08	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	410	≥ 100
Control Serum (Index)	2,7	1,9-3,5
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

1	3	4	0	3	4	7	0	0	0	0	0	0	3	8	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature 413

134034700000380000
Varicella IgM
CALIBRATOR 0.175 mL IVD
LOT 034
 2025-08 2°C-8°C
DIESSE Monteriggioni Italy CE
B2481

CONTROL SERUM

1	3	4	0	3	4	7	1	9	3	5	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature 413

134034719350000001
Varicella IgM
CONTROL 0.425 mL IVD
LOT 034
 2025-08 2°C-8°C
DIESSE Monteriggioni Italy CE
C2481

This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents 413

DATE 20.06.24

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.