

DIESE DIAGNOSTICA SENESE SpA  
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39  
Tel. +39 0577319554



### QUALITY CONTROL CERTIFICATE

KIT NAME 88000-H1 HOUSE DUST GREER

Code	88000-H1		
Lot No.	042	2025-06	2024-06

### REAGENTS

Label	Description	Lot.	Exp.	NO
D0091	88000-H1 HOUSE DUST GREER	042	2025-06	2X6
B0090	CALIBRATORE 88000-H1 HOUSE DUST GREER	042	2025-08	1
C0090	CONTROL SERUM 88000-H1 HOUSE DUST GREER	042	2025-08	1

### CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	220	≥ 100
Control Serum (KU/L)	11,4	5,5-15,6

### CALIBRATOR

6	3	5	0	4	2	6	5	1	9	4	9	3	0	8	2	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FB*

635042651949308200  
House dust-Greer Labs  
Inc.  
CALIBRATOR 0.280 mL IVD  
LOT 042  
2025-08  
DIESE Monteriggioni Italy  
CE

### CONTROL SERUM

6	3	5	0	4	2	6	2	3	4	4	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *FB*

635042623440000001  
House dust-Greer Labs  
Inc.  
CONTROL 0.450 mL IVD  
LOT 042  
2025-08  
DIESE Monteriggioni Italy  
CE

*This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Abene*

DATE *24/06/2024*

### WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.