

DIESSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554

DIESSE

Pag 1/1

QUALITY CONTROL CERTIFICATE

KIT NAME Chorus Chlamydia trachomatis IgG

Code 81261

Lot No. 042



2025-06



2024-06

REAGENTS

Label	Description	Lot	Exp.	NO
D9640	Chorus Chlamydia trachomatis IgG	042	2025-06	6X6
B9640	CALIBRATOR Chlamydia trachomatis IgG	042	2025-06	1
C9640	POSITIVE CONTROL Chlamydia trachomatis IgG	042	2025-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	537	≥ 100
Positive Control (Index)	1,5	1,1-2,0

CALIBRATOR

2 3 3 0 4 2 6 0 0 0 0 0 6 5 0 0 0 0


Label checked with barcode reader

YES



Operator Signature

[Handwritten Signature]



CONTROL SERUM

2 3 3 0 4 2 6 1 1 2 0 0 0 0 0 0 0 0 1

Label checked with barcode reader

YES



Operator Signature

[Handwritten Signature]



*This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.*

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

[Handwritten Signature]

DATE 27-06-2024

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.