


iRapid SARS-CoV-2 Quant “Neutralizing” Ab

DETERMINATION OF IMMUNE STATUS AFTER INFECTION
OR VACCINATION



iRapid SARS-CoV-2 Quant “Neutralizing” Ab

Immunochromatographic test for the **semi-quantitative** determination of IgG antibodies anti-S1 (RBD) SARS-CoV-2 virus in human serum, plasma, **fingerstick whole blood**. The test is **for professional use only** and is intended to be used as an aid in the **determination of the immune status** after infection or vaccination.

Rapid Test benefits



Fast and easy to use

Get results in a few minutes. No need to transport the biological sample to a laboratory. Easy to read results. No additional equipment is required.



Reliable

Proven diagnostics laboratory quality available on lateral flow device. Results expressed in BAU/ml standardized against 1st WHO 20/136 for anti-SARS-CoV-2. Able to detect UK, Brazilian and South African variants.



Proven Quality

Sensitivity: 96.1% CI_{95%}: 92.7 – 97.9
Specificity: 98.1% CI_{95%}: 96.5 – 98.9



Patient Comfort

Painless finger prick blood collection.

Each kit contains

- 20 cassettes individually packaged in aluminium pouches with desiccant
- 20 calibrated micropipettes to deliver up to 20 µL blood
- 20 lancets (for finger prick)
- 20 sterile alcohol wipes
- 1 buffer solution in a plastic dropper bottle



NAME	TEST	CODE
iRapid SARS-CoV-2 Quant “Neutralizing” Ab	20	70100

Check the product availability in your country by contacting your local DIESSE representative