



QUALITY CONTROL CERTIFICATE

KIT NAME MEASLES IgM CHORUS

Code	81073		
Lot No.	917	2024-06	2023-03

REAGENTS

Label	Description	Lot.	Exp.	NO
D2442	MEASLES IgM	917	2024-06	6x6
B2442	CALIBRATORE MEASLES IgM	917	2024-06	2
C2442	CONTROL SERUM MEASLES IgM	917	2024-06	2

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	470	≥ 100
Control Serum (Index)	1,9	1,4-2,6
Negative Control (Index)	0,0	< 0,7



CALIBRATOR


1	3	3	9	1	7	7	0	0	0	0	0	5	0	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *MB*

CONTROL SERUM

1	3	3	9	1	7	7	1	4	2	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES  Operator Signature *MB*



**This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.**

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *[Signature]*

DATE 07/04/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.