

QUALITY CONTROL CERTIFICATE

KIT NAME **CARDIOLIPIN-G - CHORUS**

Code	86046		
Lot No.	943	2024-06	2023-06

REAGENTS

Label	Description	Lot.	Exp.	NO
D3311	CARDIOLIPIN-G	943	2024-06	6X6
B3311	CALIBRATORE CARDIOLIPIN-G	943	2024-06	1
C3311	CONTROL SERUM CARDIOLIPIN-G	943	2024-06	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	469	≥ 100
Control Serum (GPLU/mL)	25,8	23,2-45,0

CALIBRATOR

3	3	5	9	4	3	6	3	2	1	3	7	4	5	5	8	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *FB*

336943632137465800
Cardiolipin-G
CALIBRATOR 0.176 mL IVD
LOT 943
2024-06 2°C-8°C
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

3	3	5	9	4	3	6	1	6	2	9	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature *FB*

336943616290000001
Cardiolipin-G
CONTROL+ 0.426 mL IVD
LOT 943
2024-06 2°C-8°C
DIESSSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents..... *Diavel*

DATE *11/07/2023*

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!
This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.