



DIESSSE DIAGNOSTICA SENESE SpA
53035 Monteriggioni (SI), Italy, Strada dei Laghi 39
Tel. +39 0577319554 FAX +39 0577319020

Sede Legale
20144 Milano (Italy) Via A. Solari 19 sc. 6

QUALITY CONTROL CERTIFICATE

KIT NAME CLOSTRIDIUM DIFFICILE A/B TOXINS - CHORUS

Code	81170		
Lot No.	950	2024-06	2023-07

REAGENTS

Label	Description	Lot.	Exp.	NO
D1890	CLOSTRIDIUM DIFFICILE A/B TOXINS	950	2024-06	2X6
B1890	CALIBRATORE C.DIFFICILE A/B TOXINS	950	2024-07	1
C1890	CONTROL SERUM C.DIFFICILE A/B TOXINS	950	2024-07	1
G1880	C.DIFFICILE Stool Diluent/Negative Control	C0312701	2024-07	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	360	≥ 100
Control + (Index)	3,7	2,5-4,6
Stool Diluent/Negative Control (Index)	0,2	<0,5

CALIBRATOR

9	3	6	9	5	0	6	0	0	0	0	0	2	8	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

93696060000280000
Clostridium difficile
A/B Toxins
CALIBRATOR 0.300 mL IVD
LOT 950
2024-07 2°C-8°C
DIESSSE Monteriggioni Italy CE

CONTROL SERUM

9	3	6	9	5	0	6	2	5	4	6	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

93696062646000001
Clostridium difficile
A/B Toxins
CONTROL+ 0.450 mL IVD
LOT 980
2024-07 2°C-8°C
DIESSSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements.
Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 21/07/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.