

QUALITY CONTROL CERTIFICATE

KIT NAME VARICELLA IgG - CHORUS

Code	81078		
Lot No.	947	2024-06	2023-06

REAGENTS

Label	Description	Lot.	Exp.	NO
D2471	VARICELLA IgG	947	2024-06	6x6
B2472	CALIBRATORE VARICELLA IgG	947	2024-06	1
C2471	CONTROL SERUM VARICELLA IgG	947	2024-06	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	320	≥ 100
Control Serum (Index)	4,1	2,4-4,5
Negative Control (Index)	0,0	< 0,7

CALIBRATOR

0	1	4	9	4	7	6	0	0	0	0	0	2	9	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

01494760000290000
Varicella IgG
CALIBRATOR 0.425 mL IVD
LOT 947
2024-06 2°C-8°C
DIESSE Monteriggioni Italy CE

CONTROL SERUM

0	1	4	9	4	7	6	2	4	4	5	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

01494762446000001
Varicella IgG
CONTROL+ 0.425 mL IVD
LOT 947
2024-06 2°C-8°C
DIESSE Monteriggioni Italy CE

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents.....

DATE 08/08/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.