



QUALITY CONTROL CERTIFICATE

KIT NAME TREPONEMA IgG - CHORUS

Code	81050/12		
Lot No.	952	2024-07	2023-07

REAGENTS

Label	Description	Lot	Exp.	NO
D2321	TREPONEMA IgG	952	2024-07	2X6
B2322	CALIBRATORE TREPONEMA IgG	952	2024-07	1
C2321	CONTROL SERUM TREPONEMA IgG	952	2024-07	1

CONTROLS

	Found results	Expected results
Calibrator (mOD)	370	≥ 100
Control Serum (INDEX)	2,5	2,0-3,6
Negative Control (INDEX)	0,0	<0,7

CALIBRATOR

0	1	2	9	5	2	6	0	0	0	0	0	3	6	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

01296260000360000
Treponema IgG
CALIBRATOR 0.425 mL IVD
LOT 952
2024-07 2°C-8°C
DIESE Monteriggioni Italy CE B2322

CONTROL SERUM

0	1	2	9	5	2	6	2	0	3	6	0	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

01296262036000001
Treponema IgG
CONTROL 0.425 mL IVD
LOT 952
2024-07 2°C-8°C
DIESE Monteriggioni Italy CE C2321

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents:

DATE 25/08/2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.