

QUALITY CONTROL CERTIFICATE

KIT NAME TOXOPLASMA IgA-CHORUS

Code	81043		
Lot No.	955	2024-07	2023-08

REAGENTS

Label	Description	Lot.	Exp.	NO
D2291	TOXOPLASMA IgA	955	2024-07	6X6
B2291	CALIBRATORE TOXOPLASMA IgA	955	2024-07	1
C2291	CONTROL SERUM TOXOPLASMA IgA	955	2024-07	1

CONTROLS

	Found results	Expected results
--	---------------	------------------

Calibrator (mOD)	430	≥ 100
Control Serum (INDEX)	3,2	2,0-3,7
Negative Control (INDEX)	0,1	<0,7

CALIBRATOR

0	2	8	9	5	5	6	0	0	0	0	0	3	5	0	0	0	0
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

02895660000350000
Toxoplasma IgA
CALIBRATOR 0.175 mL IVD
LOT 955
2024-07 2°C - 8°C
DIESSA Monteriggioni Italy CE 0123

CONTROL SERUM

0	2	8	9	5	5	6	2	0	3	7	0	0	0	0	0	0	1
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

Label checked with barcode reader YES Operator Signature

028956620370000001
Toxoplasma IgA
CONTROL+ 0.425 mL IVD
LOT 955
2024-07 2°C - 8°C
DIESSA Monteriggioni Italy CE 0123

This batch has been manufactured and checked according to quality requirements. Hereby, it is released to distribution.

APPROVED RESPONSIBLE HEAD, Quality Control Reagents

DATE 06-09-2023

WARNING: POTENTIAL BIOHAZARDOUS MATERIAL!

This kit may contain some reagents made with human serum or plasma. All serum or plasma used has been tested and found non-reactive for HIV-1/2, HCV and HBsAg. Because no method can offer complete assurance that HIV-1/2, HCV, HBsAg or other infectious agents are absent, reagents should be handled with maximum attention.